

# MATERIAL SAFETY DATA SHEET

NO. 129-RD-MSDS

1 OF 2

MANUFACTURER: Healthpoint Medical by DPT Laboratories

ADDRESS: 307 E. Josephine  
San Antonio, Tx 78215

PHONE: (210) 223-3281

EMERGENCY #: 1-800-441-8227

## SECTION I. MATERIAL IDENTIFICATION

A. MATERIAL NAME: Curasol Gel Wound Dressing, Hydrogel Saturated 4X4 Dressing  
TYPE: Hydrogel Wound Dressing, 4X4 Pad (gamma sterilized)

## SECTION II. INGREDIENTS

purified water, glyceryl polyacrylate, glycerin

## SECTION III. PHYSICAL DATA (APPEARANCE)

A. PRESSURIZED CONTAINER: YES  NO   
B. TYPE CONTAINER(S): LAMINATE POUCH  GLASS BOTTLE   
BOX  METAL CAN  OTHER:  
C. DESCRIPTION: CLEAR LIQUID  LIQUID EMULSION   
FRAGRANCED  CLOUDY LIQUID  CREAM   
OTHER: clear gel impregnated 4X4 pad ODOR: Faint Characteristic odor

## SECTION IV. FIRE AND EXPLOSION DATA

A. FLAMMABLE: YES  NO   
B. FLASH POINT: LESS THAN 100°F  MORE THAN 100°F   
C. CONTENTS UNDER PRESSURE: YES  NO   
D. EXTINGUISHING MEDIA: WATER  ABC ALL PURPOSE EXTINGUISHER   
CO2  OTHER:  
E. SPECIAL FIRE FIGHTING PROCEDURES: None  
F. EXPLOSIVE: YES  NO

## SECTION V. REACTIVITY DATA

A. STABLE: YES  NO   
B. HAZARDOUS DECOMPOSITION PRODUCTS: None  
C. INCOMPATIBLE WITH: None  
D. SPECIAL CONDITIONS TO AVOID: None  
E. CORROSIVE: YES  NO

## SECTION VI. HEALTH HAZARDS INFORMATION

This product is considered safe for its intended use, however, as is the case with all consumer products, some individuals may experience reactions such as redness, rash and/or swelling. Use of the product should be discontinued any time a reaction occurs.

Precautions: Avoid prolonged eye contact and oral ingestion.

First Aid:

Eyes - Remove pad, flush any adhering gel with water for 15 minutes. If irritation persists seek medical attention.

Skin - Should skin irritation, allergic reaction or rash occur, physically remove as much of the product as possible and thoroughly flush the area with water. If irritation persists, seek medical attention.

Ingestion - Remove absorbant pad from oral cavity or pharynx if possible, seek medical attention.

---

**SECTION VII. SPILL, LEAK AND DISPOSAL PROCEDURES**

---

Small Spills - Rinse the area with water.

Large Spills - Collect in a suitable container using a scoop if necessary. Material should then be disposed of by placing in a landfill or by incineration in accordance with relevant Federal, State and Local regulations.

---

**SECTION VIII. SPECIAL PROTECTION INFORMATION**

---

NONE REQUIRED

---

**SECTION IX. PRECAUTIONS AND COMMENTS**

---

A. STORAGE & HANDLING INFORMATION: Store in dry area at normal temperatures.

Avoid storing at high temperature or near heat source. Avoid freezing temperatures.

B. SPECIAL STORAGE REQUIRED: YES  NO

C. DOT CLASS: EXEMPT  RED LABEL  CONSUMER COMMODITY ORM-D   
OXIDIZER

The product described in this Material Safety Data Sheet is regulated under the Federal Food, Drug and Cosmetics Act and is safe to use as per directions on the container. The information in this data sheet has been assembled by the manufacturer based on its own studies and on the work of others. The manufacturer makes no warranties, express or implied, as to the accuracy, completeness or adequacy of the information contained herein.